

WaterPAP PRO[®]

Bubble PAP Setup
for Infants and Neonates

Accuracy

Maintains a constant
water level

Security

Depth of tube under
water held secure

Convenience

Visual representation of
prescribed pressure

Savings

Effective performance
with lower cost



www.airwaysdevelopment.com

Ordering Information on Reverse Side

Latex Free • Made in the USA • MK-19

©Copyright 2016 Airways Development LLC, All Rights Reserved
WaterPAP is a Registered Trademark of Airways Development LLC

Ordering Information

Product #	Description	Quantity
WP-7800	WaterPAP PRO Setup Includes: Coaxial Canister with 0-10cm Graduated Label, 10mm Corrugated Tube Affixed with Air Diffuser, Tube Lock Top & Protective Cover	Box of 10 Case of 40
WP-7805	WaterPAP PRO Starter Kit, Includes: Four WaterPAP Setups, One Pole Clamp, One Bracket with Tube Guide	Each
CP-2203	Pole Clamp	Each
CP-2204	Bracket	Each
WP-7704	Bracket with Tube Guide	Each
CP-2209	22mm ID x 6mm OD Fitting Combination, Connects Oxygen Supply Tube to Humidifier	Package of 40

www.airwaysdevelopment.com

Airways Development LLC
740 Jefferson Ave. Suite B
Kenilworth, NJ 07033

Phone: (908) 298-9200 • Fax: (908) 298-9201
support@airwaysdevelopment.com

